

Courtesy of : Dr. Jorge Gaspar
Instituto Nacional de Cardiología Ignacio Chávez, Mexico City, Mexico
Stent Used: Pericardium Covered Ø 4.0x18mm

Case Description:

- Globular aneurysm in the LAD preceding a 40% stenosis with a dissection flap. Stenotic plaque was dilated with a 4.0 x20 mm Rocket (ACS) balloon.

- 4.0x18mm PCS placed proximal to the take-off of the first diagonal branch with which the aneurysm was completely sealed off.

- Follow up: 6 month follow up patient continued to be asymptomatic with a < 10% luminal loss

Published: Gaspar J, Vonderwalde C, Almagor Y et al. Stent-graft (SG) with bovine pericardium (BP): acute results in porcine coronary arteries. Am J Cardiol 1998; 82:5S-6S.



Pre Intervention



Post Intervention